




Designed for Precise Balloon Dilation
with Optimal Control



KYPHOMED

Kyphoplasty Tool Kit



MEDART TECHNOLOGY PVT. LTD.

Plot No. 718, New GIDC, Gundlav, Tal & Dist. Valsad, Gujarat-396035 (INDIA)
+91 92271 53357 | info@medarttech.com | www.medarttech.com

Kyphomed System

- Performance you can count on The Kyphomed™ balloon kyphoplasty (BKP) provides MORE breakthrough solutions for you, your patients, and your procedure. Together, we can transform how vertebral compression fractures are treated.

Balloon Inflator

- This device is mainly used for hydraulic pressure on balloons to achieve the function of balloon dilation.



Features

- Optimal head fluid channel design makes it convenient for venting.
- O-type plunger design makes for low frictional resistance and can pressure easily.
- Screw fixed design of pressure gauge assures safety and reliability of the seal.
- The inflator is equipped with fluorescent gauge, which contributes to accuracy of reading.
- Ergonomic handle is comfortable and convenient to grasp.
- The inflator is of good leakproofness since it uses materials with high strength and safety.

Specifications

- Nominal Capacity 20cc
- Maximum Pressure 30atm
- Syringe Features With 10ml Syringe

Products Tool Kit



Advanced Medical-Engineering
Originating with Precision



MEDART TECHNOLOGY PVT. LTD.

Plot No. 718, New GIDC, Gundlav, Tal & Dist. Valsad, Gujarat-396035 (INDIA)
+91 92271 53357 | info@medarttech.com | www.medarttech.com